



## Mother Nature showed us how to make total fat analysis without sample transfer

To produce a perfect butterfly, Mother Nature invented the cocoon to protect the caterpillar during its transformation.

We've tried to learn by her example.

2047 SoxCap™ offers a new concept for acid hydrolysis of food and feed samples. Sample transfer in the acid hydrolysis step has, until now, been quite time-consuming for the operator. The system performs hydrolysis, filtration and washing without any sample transfer, and with minimum manual handling. Our new patented technique is based on the SoxCap capsule filter's ability to hold the fat present in the sample during hydrolysis and release it during solvent extraction in the Soxtec™ Systems.

The result is quick, simple and reliable total fat analysis – at a surprisingly low cost.

To learn more about how SoxCap sets a new standard for total fat analysis, visit [www.foss.dk/soxcap](http://www.foss.dk/soxcap) today.



**Setting new standard in total fat analysis**

Soxtec™ System 2047 SoxCap™

### Dedicated Analytical Solutions

FOSS Tecator AB  
Box 70,  
SE-263 21 Höganäs, Sweden  
Phone +46 42 36 15 00  
Fax +46 42 34 03 49  
E-mail: [info@foss.tecator.se](mailto:info@foss.tecator.se)  
[www.foss.dk](http://www.foss.dk)

Representatives in 100 countries.

**FOSS**